

**Large Area Avalanche Photodiode Detectors for a Small Etendue Thomson Scattering Diagnostic.\***

T. Biewer, M.R. Stoneking, and D.J. Den Hartog, University of Wisconsin-Madison. --- The RFP requires a close conducting shell and small field errors. For this reason diagnostic ports are kept small, leading to a small etendue (0.19 mm<sup>2</sup>-str) for TS light collection optics. In addition, the MST RFP operates at a low plasma density of  $\sim 10^{13}$  cm<sup>-3</sup>. Thomson scattered light levels are therefore lower than those found typically in tokamaks. Recently the MST TS system was upgraded to overcome the low signal level associated with its operational environment and to produce accurate single shot temperature measurements. Signal to noise enhancement was accomplished primarily through the use of 11 large area avalanche photodiode detectors with high quantum efficiency (90%) and low noise. Spectral coverage was increased from 50 nm to 100 nm. Single channel signal to noise improvement of 5X is expected, reducing the uncertainty in the electron temperature measurements from  $\sim 40\%$  to  $\sim 10\%$ . Preliminary results are reported.

\*This work was supported by the U.S. D.O.E.